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## Development of FUJI DRI-CHEM $\alpha$ -AMY Slide for the Determination of $\alpha$ -Amylase Activity in Wheat Grains

Kazuya KAWASAKI\* and Kaoru TERASHIMA\*

### Abstract

The  $\alpha$ -amylase activity in wheat grains, one of the important quality indices of wheat, can be determined with an automatic analyzer that is based on wet chemistry and that has been introduced in the field of wheat grains assessment. However, this analyzer requires time-consuming sample preparation prior to the analysis.

To provide an easier and rapider analyzing method, Fuji Photo Film has developed FUJI DRI-CHEM (FDC)  $\alpha$ -AMY Slide for the determination of the  $\alpha$ -amylase activity in wheat grains. Common to all the FUJI FDC Slides, this Slide incorporates all the reagents required for the analysis based on dry chemistry.

In the  $\alpha$ -amylase activity determination, it is desirable to conduct reaction at the optimum pH for  $\alpha$ -amylase where its activity is the maximum. But p-nitrophenol (PNP) as the detector dye exhibits a small molar extinction coefficient at that pH. To raise the sensitivity, the following two techniques have been introduced to the  $\alpha$ -AMY recipes.

(1) Reduction of the apparent pKa of PNP by using a mordant, and

(2) Optimization of the NaCl and CaCl<sub>2</sub> concentrations for the stabilization of  $\alpha$ -amylase in the extract.

The  $\alpha$ -amylase activity determination system comprising  $\alpha$ -AMY Slide described above, the FDC analyzer and related instruments were applied to some wheat samples in Hokkaido. The correlation between the data obtained with this dry chemistry system and those with the conventional Ceralpha method based on wet chemistry proved fairly good. Our dry chemistry system has various advantages including little need for complicated pre- and post-analysis operations and no waste water formation.

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
ethylidene-G7-PNP	: 4,6-Ethylidene-4-nitrophenyl-	-D-maltoheptaoside
ethylidene-G5	: 4,6-Ethylidene-	-D-maltopentaoside
ethylidene-G4	: 4,6-Ethylidene-	-D-maltotetraoside
ethylidene-G3	: 4,6-Ethylidene-	-D-maltotrioside
G2-PNP	: 4-Nitrophenyl-	-D-maltoside
G3-PNP	: 4-Nitrophenyl-	-D-maltotrioside
G4-PNP	: 4-Nitrophenyl-	-D-maltotetraoside
PNP	: 4-Nitrophenol	

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
Fig. 1 pH activity curve for malted wheat.


Fig. 2 Effect of pH on the ionization of PNP.

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
Fig. 3 Effect of the mordant on the ionization of PNP.  
Ionization rate of PNP with mordant  
Ionization rate of PNP without mordant

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
Fig. 5 Effect of NaCl and CaCl<sub>2</sub> concentrations in the extract on  $\alpha$ -amylase activity.


Fig. 6 Effect of NaCl and CaCl<sub>2</sub> in the extract on  $\alpha$ -amylase stability.

Temporal change of  $\alpha$ -amylase activity without NaCl and CaCl<sub>2</sub>

Temporal change of  $\alpha$ -amylase activity with 0.86mM $\cdot$ NaCl and 18  $\mu$ M $\cdot$ CaCl<sub>2</sub>

Temporal change of  $\alpha$ -amylase activity with 0.09mM $\cdot$ NaCl and 2  $\mu$ M $\cdot$ CaCl<sub>2</sub>

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
Fig. 7 Structure and reaction scheme of FUJI DRI-CHEM  $\alpha$ -AMY Slide.

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